

eXtreme • Power CIC  
**Destiny™**



The World's Most Powerful CIC



# The World's Most Powerful CIC




Introducing **Destiny eXtreme Power CIC**. This newest member of the Destiny family brings all of the power you'd expect from a BTE in a small package. Now you can feel confident providing your patients with severe hearing loss a discreet hearing instrument that will perform unlike any other.



## The most gain and best performing CIC available

Destiny hearing instruments with nFusion™ Technology have brought unprecedented performance to patients around the world. And now Starkey's new proprietary BluWave™ Signal Processing is driving Destiny to levels never before attainable.



*Destiny eXtreme Power CIC fitting range.*

## The industry's most advanced features in the smallest package.

### Integrated Real Ear Measurement\*

- Performs real ear measurement and incorporates data immediately throughout Inspire® OS 3.0

### Smart Features\*

- Voice Indicators – use speech to alert wearers of hearing aid status
- Self Check – runs a diagnostic check of the microphone, circuit, and receiver
- Reminder – provides an audible reminder for follow-up appointments

### Auto Path

- Auto fitting routine providing the industry's most accurate first fit

### Environmental Adaptation

- Includes Acoustic Signature, a state-of-the-art environmental detection and adaptation system
- Identifies unique sound environments and adjusts to them instantly
- Automatically detects and optimizes settings for listening in: speech, speech in noise, mechanical sounds, wind and quiet

### Data Logging

- Links with Acoustic Signature to identify and log individual sound environments
- Enables the highest level of personalization in adjustments

### Fine-Tuning

- 8 channels for sound adaptation
- 12 bands for the most precise level of fine-tuning currently available

### Active Feedback Intercept

- Virtually eliminates feedback
- Prevents whistling
- Eliminates irritating artifacts

### Automatic Telephone Response

- Instantaneously and automatically adjusts for optimum telephone communication

### Adaptive Indicator Tones

- Unique tones for memory, low battery, volume, etc.
- Automatically adjusts in noisy environments

\*Available only on the Destiny 1600 series.



A global hearing technology company  
headquartered in Eden Prairie, Minnesota, U.S.A.

Starkey Laboratories, Inc.  
6700 Washington Avenue South  
Eden Prairie, MN 55344-3476  
800.328.8602

[StarkeyPro.com](http://StarkeyPro.com)