

AS110

Preliminary Screening Device

AUDITORY SCREENER



- LIGHTWEIGHT
- DURABLE
- ECONOMICAL



AS 110 AUDITORY SCREENER

The lightweight and durable AS 110 is an excellent preliminary screening device. The AS 110 features pure tone air conduction testing and offers increased control and flexibility. Weighing only nine lbs. and housed in an impact resistant case, the AS 110 is ideal for individual or group testing in schools, doctors' offices and industrial environments. This portable unit comes standard with TD-39 headphones.

Additional features include:

- Pure Tone Test Frequencies of 250-8000 Hz
- Hearing Threshold Dial from -10 to 110 dB HL in 5 dB intervals
- Continuous, Pulsed or Warble Tones
- Reversible Tone Interrupter Switch
- Patient Signal
- Indicator Lights for Earphone Selection and Pure Tone Presentation

Technical Specifications AS 110 Auditory Screener Audiometer

Power Source: 110 or 220 volts A.C. only; 50-60 Hz

Power Consumption: 19 watts

Frequency Range (Air): 125-8000 Hz

Accuracy of Frequency Calibration: Less than 3% error from indicated frequency

Distortion: Meets or exceeds ANSI S3.6 - 1996

Tone Interruption: Per ANSI S3.6 - 1996

Accuracy of Intensity: Signal level at any dial setting is +/- 3 dB of reference SPL from 125 to 4000 Hz, and +/- 5 dB at 6 K and 8 K Hz

Weight in Pounds: 9

Portable Case Dimensions: 13" wide x 13" deep x 5 3/4"

Group Presentation Capability: At all frequencies with up to 30 phones

Additional Features: Pulsed or continuous tone presentation, reversible tone presentation, function switch indicator, lights, patient signal, warble tone (FM)

STARKEY'S IRON-CLAD WARRANTY

Each STARKEY AUDIOMETER is warranted against defects in workmanship and material as to parts arising in connection with normal usage for a period of one (1) year from date of purchase.

STARKEY'S obligation under the above warranty is limited to replacing or repairing any defective part or parts of the audiometer for the one (1) year duration of the warranty.

